




www.DelpMixer.com
Info@DelpMixer.com
1-800-504-1935 Toll Free

2528 Billingsley Road
Columbus, Ohio 43235
614-889-7331 (F)

OVA10 Open Vessel Adapter for Standard 10mm Shafting



Open Vessel Adapter 10 Assembly

Includes:

Flexible Drive Cable, 36 or 60 inch length (attaches to any drive motor with a 5/16 (8mm) diameter shaft), OVA 10 Stand Mount and OVA10 Adapter

- Non-sparking device that eliminates explosion hazard when mixing volatile materials
- Suitable for use in all open vessels or beakers
- Operates with all chuck style mixers, including IKA "Eurostar," Troemner "TallBoy," or "Stir-Paks"
- Accepts standard 10MM glass, coated or metal shafts and standard agitators
- Manufactured from high-grade materials to withstand autoclave sterilization